

Table S1. Comparison of web-based oligonucleotide design tools. The comparison between different oligonucleotide design/check tools shows that most of the tools did not involve specificity and cross dimer checks. Therefore, these tools also don't involve the possibility for multiplexing. Except for oli2go, which includes all listed design principles, MFEprimer involves most of the processing steps. However, it only provides limited access to background data for specificity checks and less comprehensive cross dimerization calculations.


Tool	Primer Design	Probe Design	Hairpin Check	Batch design*	Probe Specificity Check	Primer Specificity Check	Cross Dimer Check	Multiplexing
Oli2go	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
MSP-HTPrimer	Yes	Yes	Yes	Yes	No	No	No	No
PrecisePrimer	Yes	No	No	Yes	No	No	No	No
Primer3	Yes	Yes	Yes	No	No	No	No	No
MFEprimer	No	No	Yes	Yes	Yes (limited)	Yes	Yes	Yes (limited)
Primer-BLAST	Yes	No	No	No	No	Yes	No	No

* Batch design means the application of the software on more than one input sequences in one run. However, this term does not cover multiplexing, as design steps such as specificity and secondary structure checks need to be performed for multiplex applications.

Table S2. Comparison between cross dimer checks using oli2go and MFEprimer. Primer sequences of 6 genes (ceoA: fwd – TTCAAGGACGGCGCG, rev – TATAGCCGAGRTTGATGCG; cphA4: fwd – GCGAGCTGCACAAGCT, rev – CCTTGCGGGTAAAGGC; cfiA7: fwd – TTATCCTTATCTCCATGCTT, rev – TTCGGCGAGGGATACATAAGT, OXY2-5: fwd – GTGCAGCACCAGTAAAG, rev – CGTTAATCTCCAGCCTTT; ykkC: fwd – TTAACATGGAGCGGCACT, rev – TTATGCCTCGCCTCCT; ykkD: fwd – ATGCTGCACTGGATCAGTTTA, rev – GCAAACCAACAATGATCAACAG) were used to compare cross dimerization check results of MFEprimer and oli2go. Parameters for both tools were chosen as follows: concentration of monovalent cations 10 mM, concentration of divalent cations 22 mM, concentration of dNTPs 1.75 mM, primer concentration 50 nM, deltaG threshold -10.000 cal/mol and a T_m threshold of 40 °C. First, MFEprimer does not involve the possibility to define thresholds for deltaG and the T_m of secondary structures. Therefore, resulting secondary structures must be examined manually. If we apply the same threshold for MFEprimer as used for oli2go, there would be no cross dimerization resulting from MFEprimer. However, oli2go would detect 2 primer dimers using the same parameters as with MFEprimer. This difference in results of both software packages arises from the fact, that MFEprimer involves less thermodynamical parameters and methods for finding the most stable structure. Furthermore, oli2go compares the melting temperature and the delta G value of the secondary structure. If only one of these parameters exceeds the threshold, the cross dimer check fails. However, MFEprimer only returns delta G values of the secondary structures. These characteristics lead to a less stringent evaluation of secondary structures.

Tool	Primer 1	Primer 2	DeltaG	T _m	Dimerisation
Oli2go	cfiA7 rev primer	ykkC ref primer	-11324.5	37.43	True
	OXY2-5 fwd primer	ykkD fwd primer	-10246.7	32.51	True
MFEprimer	cphA4 rev primer	OXY2-5 rev primer	-7600.0	-	-

Figure S1. Agarose gel image of experimental primer evaluation experiments. Primers were evaluated in single-plex PCRs to analyze the functionality and specificity of the respective primer pairs. The PCR was conducted using the DNA-free Mastermix 16S Basic PCR kit (Molzylm, Germany) according to the manufacturer's instructions comprising 200 nM of each primer, 1 µl of target DNA and 0.32 µl of the Moltag 16S polymerase in a total volume of 10 µl. The thermal cycling was performed as follows: 94°C for 5 min; 40 cycles of 94°C for 30 s, 46°C for 30 s, 72°C for 30 s; and a final elongation cycle at 72°C for 7 min. The highlighted numbers indicate deviations from the expected amplification product size or low amplification efficiencies. In such cases, the primers have to be redesigned if experimental mistakes can be excluded.



REFERENCES

1. Barišić, I., Kamleithner, V., Schönthaler, S. and Wiesinger-Mayr, H. (2014) Fast and highly specific DNA-based multiplex detection on a solid support, *Applied microbiology and biotechnology*.